

OPTOCAM-I

5MPx CMOS Microscope Camera - SpecSheet

- USB CMOS Color Camera
- Capturing of photos and videos
- 1/2.5 " CMOS sensor, progressive scan
- C-Mount connection
- Resolution up to 2592 x 1944 Pixel (6 FPS)
- Up to 52 images/sec @ 640x480
- Dynamic range 8 Bit
- Exposure times from 1/10000s up to 30s
- USB 2.0 (cable included)
- Capturing software "OptoCapture" included
- Supported OS:
 - Windows XP-SP3 (32bit)
 - Windows Vista (32bit)
 - Windows 7 (32bit/64bit)
- Software languages: english, german
- USB-stick with software & manual included

